




High-throughput behavioral testing systems that enable real-time measurements of neural activity with unprecedented precision


Drinko MEASURER

For accurate Fluid Intake Measurements



Features

- ▶ Dual ball bearings inside the tube restrict water output and prevent excess leakage.
- ▶ Choose between tube and bottle lengths.
- ▶ All parts are autoclavable for easy cleaning.
- ▶ Ideal for conditioned taste aversion, drug dosing and alcohol consumption experiments.



Tube Length	65 mm (2.5 in) 90 mm (3.5 in) 120 mm (4.5 in)
Clip	Pen type clip
Material	Body of cap: silicon Tube: stainless (SUS304) Pen type clip (SUS304)



Volume	15 mL bottle 10 mL bottle
Material	Heat-resistant glass

OB04-0919

info@amuzainc.com | (858) 225-6869